

BRAIN-PAD "DOUBLE" MOUTHGUARD

DO NOT USE this mouth guard IF you have an under-bite (Class III prognathic jaw relationship).



The BRAIN-PAD "DOUBLE" is a LATEX-FREE product!

Custom Fitting Instructions:

FITS OVER BRACES: Immediately After boiling, channel flanges may be widened with thumb. After creating impressions, DO NOT push outside flanges inward against braces/teeth.

OPTIONAL SIZING:

STANDARD MODEL
Do not cut in front of line C.

JUNIOR MODEL
Do not cut in front of line B.

If unit is cut for re-sizing -

Time critical in seconds.

Time critical in seconds.

2 STAGE BOILING PROCESS

Completely insert Boil-N-Bite Stick as pictured:

FIRST STAGE:
20 SECONDS

Place in water up to boil line.

DO NOT OVER BOIL



SECOND STAGE:
10 SECONDS

Place in water completely.

OVER 30 SECONDS TOTAL DESTROYS THE UNIT!



AFTER BOILING:

- 1) LEAVING BITE STICK IN BREATHE SPACE, PLACE CLENCHER ON LOWER TEETH FIRST.
- 2) BRING UPPER TEETH INTO TOP CHANNELS.
- 3) BITE FIRMLY FOR 10 seconds.
- 4) RUN UNDER COLD WATER TO SET IMPRESSIONS. REMOVE BITE STICK.

Sanitize regularly with an astringent mouthwash for good hygiene. With proper care, your guard should last for a whole season. 1.

NON-STRAP Models:
(DG/DGY-200)
DOES NOT HAVE strap attachment Tab.

STRAP Models
(DG/DGY-205)
DOES HAVE Strap Attachment Tab built into front of unit.

STRAP MODEL
DG/DGY-205

CORRECT STRAP ATTACHMENT
Inside/behind facemask bar.

WARNING:
DO NOT Knot, Tie, or Wrap any part of string attachment.
Doing so can negate Safety Clasp's functionality & effectiveness!

[detach along perforation and KEEP below information to satisfy warranty requirements]

WARRANTY by Brain-Pad, Inc. - (Hereinafter referred to as "the company")

The company's liability for a defective 'DOUBLE' mouthguard(DG) is expressly & solely limited to (i) the replacement of the DG; or (ii) the lesser of (a) actual dental expenses incurred by the Purchaser; or (b) the sum of Two Hundred Three (\$203) per injured or replaced tooth, up to a maximum liability of Six Thousand Five-hundred Dollars (\$6,500) per incident; provided, however, the company's obligations under the immediately preceding sentence shall be reduced by the amount of any insurance collected is payable to Purchaser for the damage, loss, and/or injury claimed by Purchaser on account of the DG. Notwithstanding any other provision in the Warranty to the contrary, (i) this Warranty does not provide or shall be deemed to provide payment or indemnity for orthodontic treatments or preparation of the tooth to be used as an anchor for a bridge, (ii) the company shall have no responsibility or liability whatsoever under this Warranty unless (a) the DG was fitted by Purchaser in accordance with the printed instructions accompanying this Warranty; and (b) the original Instructions/Warranty Document received with the purchased DG and a complete written notice from Purchaser of the injury, loss and/or damage is received by the company within thirty (30) days after the occurrence of the injury, loss, and/or damage ("Notice") addressed to Brain-Pad, Inc., Warranty Department, 322 Fayette St., Conshohocken, Pennsylvania, 19428.

Other limitations apply. View complete warranty information at: www.brainpads.com.

Manufactured exclusively by Brain-Pad, Inc., 322 Fayette St., Conshohocken, PA 19428.
610.397.0893 www.brainpads.com Copyright 2000-2007 Brain-Pad, Inc.
Protected under U.S. & Foreign Patents:

5,636,379., # D493,578S., # D500,895S., # D523,994., # D525,749.,
D526,093., # D526,095., # D527,848., # D530863., # D532,559., #D537,986S.,
#D537,987S., # 000468947-0001., # 000468947-0002., # CA2251672.

All Rights Reserved. Other U.S. & Foreign Patents pending. Made in Taiwan.

ADULT WITH STRAP



7 57920 90020 9

DG-205 Black

INST-08 07